

How to isolate nondenatured protein extracts from Neurospora for western analysis.

Yi Liu

Procedure:

Reagents

Extraction Buffer: 50mM HEPES (pH7.4), 137 mM NaCl, 10% Glycerol.

Add protease inhibitors fresh to the buffer just before the extraction.

Protease inhibitors final concentration:

Pepstatin A (1mg/ml), Leupeptin (1mg/ml), PMSF (1mM)

Stock kept in -20°C:

Pepstatin A: 1mg/ml in methanol.

Leupeptin: 1mg/ml in H₂O.

PMSF: 100 mM in isopropanol.

For certain experiments with PPase inhibitors (made fresh):

25mM NaF, 10mM Na₄P₂O₇.10H₂O, 2mM Na₃VO₄, 1mM EDTA.

1. Grind tissue in liquid nitrogen with a mortar and pestle.
2. Suspend ground tissue in ice-cold extraction buffer containing protease inhibitors. Mix by gentle inversion. Let sit on ice for 5 min and then mix again.
3. Centrifuge cell homogenate at 12,000 rpm for 15 min at 4°C.
4. Keep the supernatant, and measure the protein concentration.
2 ml protein extract to 800 ml H₂O, then 200 ml Bio-Rad protein assay dye. Mix and measure OD₅₉₅ and calculate protein concentration.
5. Prepare protein sample for SDS-PAGE. For small gel (Bio-Rad Protean III gel) ~20-40mg per well.

	Small gel (50mg)
H ₂ O	(20-x) ml
2X loading buffer	20 ml
protein extract	x ml

6. Boil the sample on heat block for 5min, chill on ice and load onto SDS-acrylamide gel.

Urea containing 2X loading buffer (the addition of Urea is optional):

100mM Tris pH6.8, 4% SDS, 2mM EDTA, 2% glycerol, 6M Urea.

To use: mix 900 ml buffer, 50 ml 2M DTT and 50ml 5% Bromophenol Blue.

References:

Garceau NY, Liu Y, Loros JJ, Dunlap JC. Alternative initiation of translation and time-specific phosphorylation yield multiple forms of the essential clock protein FREQUENCY. *Cell*. 1997, 89(3): 469-76.