

How to obtain naked nuclei extruded from the cell.

Namboori B. Raju

Background

Naked nuclei extruded from conidia or other cell types are potentially useful for studies in cell biology.

Procedure

The following is copied from Raju (1998):

The silver-staining procedure described by Lu (1993) for surface-spreading synaptonemal complexes of *Neurospora* is also suitable for staining expelled paraphysal nuclei. Each nucleus shows a darkly-stained prominent nucleolus. I have adapted the staining method for routine observation of conidial nuclei. Conidia were fixed in a few drops of 4% paraformaldehyde for 3 or 4 min at 4°C. A drop of conidial suspension was then placed on a glass slide and squashed under a cover glass with gentle blows from a rubber hammer to break open the conidia. The cover glass was lifted off after freezing the slide on a Cold Plate and the nuclei were flooded in a few drops of 4% paraformaldehyde + 0.04% SDS for 15 to 20 min. The slide was rinsed in distilled water, dipped in 0.2% Kodak Photoflo, and air-dried. For silver staining, one drop of 50% silver nitrate and one drop of gel developer (2% gelatin + 1% formic acid) were mixed together and added on to the slide and covered with a cover glass. The slide was heated on a hot plate at 45°C for 2-3 min during which the nuclei were stained brownish to dark brown. The cover glass was promptly floated off and the slide was rinsed in distilled water and mounted under a new cover glass in a drop of 10% glycerol for microscopic examination. I find the silver staining particularly useful for examining nucleolus behavior in the segmental duplication strain *Dp(AR33)* that contains two normal nucleolus organizers (Perkins *et al.* 1980), and in the translocation strain *T(OY321)* that contains a split nucleolus organizer (Perkins *et al.* 1984). In both strains, up to 10% of nuclei show two nucleoli (slightly smaller than normal) while the remaining nuclei contain a single larger nucleolus, interpreted as resulting from fusion. The extracted conidial nuclei may also be used for *in situ* hybridization studies.

References

- Lu, B. C. 1993. Spreading the synaptonemal complex of *Neurospora crassa*. *Chromosoma* 102: 464-472.
- Perkins, D. D., N. B. Raju, and E. G. Barry. 1980. Chromosome rearrangement in *Neurospora* that produces viable progeny containing two nucleolus organizers. *Chromosoma* 76: 255-275.
- Perkins, D. D., N. B. Raju, and E. G. Barry. 1984. A chromosome rearrangement in *Neurospora* that produces segmental aneuploid progeny containing only part of the nucleolus organizer. *Chromosoma* 89: 8-17.

Raju, N. B. 1998. A quick silver-staining method for conidial nuclei, and its usefulness for studying the behavior of the nucleolus organizer in *Neurospora*. *Neurospora 1998* (Asilomar conference), Abstract 34, p. 25.

NBR