

## How to amplify DNA from macroconidia

Susan Crosthwaite and Cas Kramer

### Background

To test the genotype of a strain or to amplify a region of genomic DNA for cloning or sequencing it is not always necessary to first extract and purify genomic DNA from other cellular constituents. The region of interest can be amplified directly from a suspension of asexual spores. *Neurospora* produces three distinct types of spores: ascospores, macroconidia and microconidia. Ascospores are the product of a sexual cross between strains of opposite mating type<sup>1</sup>. Macroconidia and microconidia are asexual spores. Macroconidia contain between 1-5 nuclei and are produced in abundance on most media. Microconidia are much less abundant, smaller than macroconidia and contain most often one, but occasionally two nuclei<sup>2</sup>.

### Procedure

For *Neurospora* template DNA harvest macroconidia (5-10 days old) with a loop or sterile toothpick from slants. Transfer to 100 µl sterile dH<sub>2</sub>O and vortex. Dilute 1µl in 999 µl dH<sub>2</sub>O and measure the OD<sub>530</sub> reading to estimate the number of spores using the following formula:

$$\text{OD}_{530} \times \text{dilution factor} \times 10^6 = \text{conidia/ml}$$

\*For amplification use  $2 \times 10^4$  spores/ 50 µl reaction. Incubate the spore suspension at 30 °C for 75 minutes then freeze in liquid nitrogen.

dH <sub>2</sub> O	31.75 µl
10 x PCR buffer	5µl
macroconidia * ( $10^6$ spores/ ml)	2µl
10 mM dNTP mix	1 µl
forward primer 50 pmol/ul	5 µl
reverse primer 50 pmol/ul	5 µl
taq polymerase (5 units/µl)	0.25 µl
Total volume	50 µl

The first denaturation step should be at least 5 minutes at 95 °C. Subsequently denature for 30 to 60 seconds at 95 °C. Carry out 40 cycles of denaturation, annealing of primers, and extension.

Henderson *et al.*<sup>3</sup> have reported a slightly different method from the one described above. A loop full of conidia is suspended in 100 µl Tris-EDTA (pH 8.0), placed in a boiling water bath for 5 minutes then transferred to ice. The sample is then placed in a

microcentrifuge and spun at 13 K rpm for 5 minutes to pellet the cell debris. 5  $\mu$ l of the supernatant is used per 50  $\mu$ l PCR. They note that successful amplification of the genomic DNA is dependant on the polymerase used. In their hands Red Hot DNA polymerase from ABgene works well.

### **References**

1. Shear, C. L. and Dodge, B. O. (1927) Life Histories and heterothallism of the red bread mold fungi of the *Monilia sitipholoia* group. J. Agric. Res. 34, 1019-1042.
2. Dodge, B. O. (1932) The non-sexual and sexual functions of microconidia of neurospora. Bull. Torrey Botan. Club. 59, 347-360.
3. Henderson S. T., Eariss G. A. and D. E. A. Catcheside (2005) Reliable PCR amplification from *Neurospora crassa* genomic DNA obtained from conidia Fungal Genet. Newsl. 52: 24.