

# Immunodepletion

Yi Liu

**This is for studying quantitative relationships of proteins in a complex.**

## Procedure

1. Couple the antisera to protein G beads:
  - a. Dilute 100  $\mu$ l of serum in 1 ml of IP buffer (recipe) with 10mg/ml of BSA and 0.1% NP40.
  - b. Add 50  $\mu$ l of protein G beads (packed volume), rotate on ? at 4 °C overnight.
  - c. Wash the beads with X amount of the IP buffer 3 times.
2. Deplete:

Mix 1mg total extract diluted in 1ml IP buffer + 20 $\mu$ l protein G beads + 0.05% NP40, rotate on what overnight at 4 °C
3. Pellet the protein G beads with a brief spin at 4 °C, 4,000g X 1min, take the supernatant and load 20  $\mu$ l onto SDS-PAGE for western blot.

Reference: Cheng et al (2005) Genes & Development 19:234-241